

## Product Information

### SIGMAFAST™ 3,3'-Diaminobenzidine tablets

tablet, To prepare 5 mL

Catalog Number **D4293**

Storage Temperature  $-20\text{ }^{\circ}\text{C}$

#### Product Description

3,3'-Diaminobenzidine is used in many applications for visualization of peroxidase activity.<sup>1-6</sup> SIGMAFAST™ 3,3'-Diaminobenzidine (DAB) tablets have been developed for use in immunohistology as a precipitating substrate for the detection of peroxidase activity. DAB is the immunohistology substrate of choice, as it produces an intense brown-black stain, which is resistant to alcohol. Slides stained with DAB may be coverslipped in the traditional manner and stored for future reference.

Peroxidase + 2 H<sub>2</sub>O<sub>2</sub> → O<sub>2</sub> + 2 H<sub>2</sub>O (pH 7.6)

O<sub>2</sub> + DAB → insoluble, brown-black precipitate

SIGMAFAST DAB tablets require no additional ingredients or procedures to prepare an active substrate solution. One SIGMAFAST DAB tablet set (one DAB tablet and one Urea Hydrogen Peroxide tablet) dissolved in 5 mL of ultrapure water provides 5 mL of ready-to-use substrate solution that contains:

3,3'-Diaminobenzidine (DAB)	0.7 mg/mL
Urea Hydrogen Peroxide (H <sub>2</sub> O <sub>2</sub> equivalence, 0.7 mg/mL)	2.0 mg/mL
Trizma® buffer	60 mM

This product has been used in studies on such systems as hepatocellular carcinoma,<sup>7</sup> gastric carcinoma,<sup>8</sup> *Xenopus laevis* muscle and embryos,<sup>9</sup> mouse models,<sup>10-14</sup> *Zootoca vivipara* embryos,<sup>15</sup> dogfish testes,<sup>16</sup> rainbow trout,<sup>17</sup> Sprague-Dawley rats,<sup>13</sup> and food samples.<sup>18</sup>

#### Components

3,3'-Diaminobenzidine (DAB) tablets Catalog Number D9292	5 or 50 set
Urea Hydrogen Peroxide tablets Catalog Number U1380	5 or 50 set

#### Reagents and Equipment Required but Not Provided

- Ultrapure water
- Pipette capable of delivering 5 mL
- Test tubes
- 0.2 μm filter (e.g., Catalog Number WHA10462701)
- Nickel(II) chloride hexahydrate (NiCl<sub>2</sub>, e.g., Catalog Number 223387) or Cobalt(II) chloride hexahydrate (CoCl<sub>2</sub>, e.g., Catalog Number 202185), to prepare 0.3% (w/v) stock solution for enhancement of tissue stains
- PBS (e.g. Catalog Number P4417) for washing

#### Precautions and Disclaimer

This product is for Research Use Only. Not for Use in Diagnostic Procedures. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

#### Preparation Instructions

1. Remove the required number of DAB and Urea Hydrogen Peroxide tablets from the freezer.
2. Allow the tablets to reach room temperature.
3. Open DAB tablet package (silver foil) and Urea Hydrogen Peroxide tablet package (gold foil) and drop the tablets into an appropriate container. **Do not touch the tablets with your fingers.**
4. Add 5 mL of ultrapure water.
5. Vortex until dissolved.

The SIGMAFAST DAB Substrate Solution is now ready for use. For best results, the solution should be used within one hour.

#### Storage/Stability

Store the tablets at  $-20\text{ }^{\circ}\text{C}$ .

### Procedure

1. Cover the treated tissue section with 0.2–0.5 mL of DAB Substrate Solution.
2. The DAB reaction may occur rapidly. Color development should be carefully monitored during the reaction to prevent overdevelopment and high backgrounds. Reactions may be stopped by gently washing the slide in water or PBS.
3. DAB reactions may be enhanced by the addition of a NiCl<sub>2</sub> or CoCl<sub>2</sub> solution. Add 0.5 mL of 0.3% (w/v) stock solution to 4.5 mL DAB Substrate Solution. The addition of metal salts to DAB changes the color of the precipitate product to black or blue-black.
4. Occasionally the DAB Substrate Solution may be hazy. The haziness may be removed by filtering the solution through a 0.2 µm filter.
5. Tissues stained with SIGMAFAST DAB Substrate Solution may be dehydrated with alcohol and mounted with traditional resinous mounting media.

### Troubleshooting

#### Background is too high

1. Use a blocking step prior to the application of the primary antibody. Diluted normal serum (10% v/v) from the same species as the secondary antibody generally produces the best results.
2. Block endogenous peroxidase by flooding the slide with a solution of 4 parts methanol and 1 part 3% H<sub>2</sub>O<sub>2</sub> solution.
3. Decrease the staining time.
4. Titer the conjugate to optimize working dilution.

#### No color develops or color is too faint

1. Adjust the concentration of the primary antibody.
2. Adjust the concentration of the secondary antibody.
3. Determine if the enzyme conjugate is active.
4. Consider using an amplifying system such as avidin-biotin or peroxidase anti-peroxidase.
5. Increase the staining time.
6. Determine if enzymatic treatment (unmasking) of the antigen is required prior to application of the primary antibody.

### References

1. Nakane, P.K., and Pierce, G.B., Jr., *J. Histochem. Cytochem.*, **14(12)**, 929-931 (1966).
2. Trojanowski, J.Q. *et al.*, *J. Histochem. Cytochem.*, **31**, 1217-1223 (1983).
3. DeJong, A.S.H. *et al.*, *Histochem. J.*, **17(10)**, 1119-1130 (1985).
4. Chu, N.M. *et al.*, *J. Histochem. Cytochem.*, **37(2)**, 257-263 (1989).
5. Merchenthaler, I. *et al.*, "Silver Intensification in Immunocytochemistry", in *Techniques in Immunocytochemistry* (Bullock, G., and Petrusz, P., eds.). Academic Press Ltd. (San Diego, CA), pp. 217-252 (1989).
6. Hsu, S., and Soban, E., *J. Histochem. Cytochem.*, **30(10)**, 1079-1082 (1982).
7. Yarchoan, M. *et al.*, *Clin. Cancer Res.*, **23(23)**, 7333-7339 (2017).
8. Tsamandas, A.C. *et al.*, *Anticancer Res.*, **29(2)**, 703-709 (2009).
9. Lunardi, A., and Dente, L., *Mech. Dev.*, **119 Suppl 1**, S49-S54 (2002).
10. Nolan-Stevaux, O. *et al.*, *Genes Dev.*, **23(1)**, 24-36 (2009).
11. Brotherton, T.E. *et al.*, *Proc. Nat. Acad. Sci. USA*, **109(14)**, 5505-5510 (2012).
12. van Dop, W.A. *et al.*, *Gastroenterology*, **136(7)**, 2195-2203 (2009).
13. Nandi, S.S. *et al.*, *Diabetes*, **65(10)**, 3075-3090 (2016).
14. Wakabayashi, T. *et al.*, *Arterioscler. Thromb. Vasc. Biol.*, **38(11)**, 2576-2589 (2018).
15. Stewart, J.R. *et al.*, *J. Exp. Biol.*, **214(18)**, 2999-3004 (2011).
16. Gautier, A. *et al.*, *Biol. Reprod.*, **91(4)**, 1-15 (2014).
17. Jørgensen, T.R. *et al.*, *Folia Parasitol. (Praha)*, **56(4)**, 251-258 (2009).
18. Boutten, B. *et al.*, *Food Agr. Immunol.*, **11(1)**, 51-59 (1999).

SIGMAFAST is a trademark of Sigma-Aldrich Co., LLC. Trizma is a registered trademark of Sigma-Aldrich Co., LLC.

PCS,RGB,GCY,MAM 11/18-1