

Research Report

Extended, High Density Growth of CHO-K1 Cells in EX-CELL® 302 Serum-Free Medium

Introduction

SAFC® has developed a serum-free medium, EX-CELL 302, which meets or exceeds requirements for use with CHO-K1 cells grown in suspension culture. Raw materials (amino acids, lipid supplements and carbohydrates) used in this medium are of the highest quality: USP, EP, JP or ACS, and are non-animal sourced wherever possible. With the addition of 4 mM L-glutamine, EX-CELL 302 will support growth of CHO-K1 cells in batch suspension cultures for greater than eight (8) days while maintaining high cell viability (> 70%). This long culture period and stationary phase allow for an extended production period and higher culture productivity. L-glutamine, hypoxanthine and thymidine have been omitted from EX-CELL 302 allowing for use with either the DHFR^r or GS selection systems.

Materials and Methods


- CHO-K1 cells (ATCC CCL 61) were obtained from the American Type Culture Collection (Rockville, MD) and were previously adapted to EX-CELL 301 (Catalog No. 14331C-Available for custom development only).
- L-Glutamine (Catalog No. 59202C), SAFC, Lenexa, Kansas.
- EX-CELL 302 (Catalog No. 14312C-Available for custom development only), SAFC, Lenexa, Kansas.
- EX-CELL 301 (Catalog No. 14311C-Available for custom development only), SAFC, Lenexa, Kansas.

Cell Culture

- Cultures were maintained in an incubator in 250 mL Bellco spinner flasks using a seeding density of 2×10^5 cells/mL in 100 mL of EX-CELL 302. Impeller speed was set to approximately 60 rpm.
- Daily cell counts and viabilities were determined by Trypan Blue dye exclusion.


Results

Figure 1. EX-CELL 302 CHO-K1 Growth Profile (10 passages)



CHO-K1 cultures were monitored daily for growth and viability. These data represent the average cell density and viability for 10 passages. Cultures were seeded at a density of 2×10^5 cells/mL in 100 mL media using a spinner speed of 60 rpm. An increase in growth rate and final cell density were seen in subsequent passages after full adaptation.

Figure 2: EX-CELL® 302 CHO-K1 Multiple Passage



Spinner cultures were initiated from stationary cell cultures that were grown in EX-CELL 301 medium. Cultures were passaged every 4 - 5 days using a seeding density of 2×10^5 cells/mL. Cultures were monitored daily for maximum cell density after passaging. Cultures were terminated when culture viability dropped below 70%.

Conclusions

Data presented here demonstrate extended cellular growth in EX-CELL 302 medium for eight (8) days in spinner batch cultures. Cultures of CHO-K1 cells achieved an average maximum cell density of approximately 3.5×10^6 cells/mL with exponential doubling times in the range of 20 hours while maintaining viabilities greater than 70%. EX-CELL 302 contains low levels ($< 100 \mu\text{g/L}$ total) of recombinant proteins. CHO-K1 cells, cryopreserved and adapted to serum-free growth in EX-CELL 301, were transferred directly into EX-CELL 302 without an adaptation period. SAFC's regulatory friendly EX-CELL 302 is suitable for the production of human therapeutic products since it does not contain animal-derived proteins or other undesirable components.

For more information, visit safcglobal.com or contact SAFC

USA: + 1-913-469-5580 • info-na@sial.com

Europe: +44 (0)1264-333311 • info-eu@sial.com

Asia/Pacific: +61 (0)2 9841 0555 • auscustserv@sial.com

Warranty, Limitation of Remedies

SAFC warrants to the purchaser for a period of one year from date of delivery that this product conforms to its specifications. Other terms and conditions of this warranty are contained in SAFC written warranty, a copy of which is available upon request. ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE IMPLIED WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE EXCLUDED. In no case will SAFC be liable for any special, incidental, or consequential damages arising out of this product or the use of this product by the customer or any third party based upon breach of warranty, breach of contract, negligence, strict tort, or any other legal theory. SAFC expressly disclaims any warranty against claims by any third party by way of infringement or the like. THIS PRODUCT IS INTENDED FOR PURPOSES DESCRIBED ONLY AND IS NOT INTENDED FOR ANY HUMAN OR THERAPEUTIC USE.

Additional Terms and Conditions are contained in the product Catalog, a copy of which is available upon request.

©2012 Sigma-Aldrich Co. LLC. All rights reserved. SAFC and SIGMA-ALDRICH are trademarks of Sigma-Aldrich Co. LLC, registered in the US and other countries. EX-CELL is a trademark of Sigma-Aldrich Co. LLC, registered in the US and other countries. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see product information on the SAFC website at www.safcglobal.com and/or on the reverse side of the invoice or packing slip.

Issued June 2012 R003
 0103 1205 0306