


cis-Diammineplatinum(II) dichlorideProduct Number **P4394**
Store at Room Temperature

CAS# 15663-27-1

Synonyms: Cisplatin; cis-DDP; CPDD; CDDP;
NCI-C55776; Neoplatin; NK 801; NSC-119875;
Peyrone's Chloride**Product Description**

Appearance: Yellow to orange powder

Molecular Formula: Pt(NH₃)₂Cl₂

Molecular Weight: 300.0

Melting Point: 270 °C (decomposes)¹

Extinction Coefficient: (0.1 N HCl):

At 203 nm, E^{mM} = 5.2 mM⁻¹cm⁻¹, E^{1%} = 173;301 nm, E^{mM} = 0.130 mM⁻¹cm⁻¹, E^{1%} = 4.33;362 nm, E^{mM} = 0.024 mM⁻¹cm⁻¹, E^{1%} = 0.806.¹

A method of preparation and physical properties including the pK_a of several species of cisplatin in water,² the UV and visible spectra, HPLC methods, and mass spectra have been reported.¹

Cisplatin is a platinum-containing, broad activity antineoplastic and alkylating agent effective against solid tumors of testes, ovaries, and bladder; epithelial malignancies; and cancers of the esophagus, lung, head, and neck.^{1,3,4} Cisplatin is effective against sarcoma and leukemia in mice.⁵ Cisplatin enters cells by diffusion. Its chloride is replaced with water forming the active, positively charged species, which may react with DNA to form intra (between N7 atoms of adjacent pairs) and interstrand crosslinks, thereby, inhibiting DNA replication.⁴ Inhibition of the BamH I cleavage and unwinding of pBR322 DNA by cisplatin was reported.⁶

Cisplatin (12.5 μM) binds to free sulfhydryl groups in tubulin and causes a partial depolymerization of microtubules.⁷ It changes microtubule assembly by

direct tubulin modification⁸ and causes changes in the cytoskeletal pattern of tumor cells.⁹ Cisplatin has induced apoptosis in cisplatin-sensitive and -resistant human ovarian cancer cells.^{10,11} Cisplatin (5 mM) complexes with phosphatidylserine lipids in the membrane¹² and reportedly induces nephrotoxicity in the rat mainly by proximal tubular impairment.¹³

The mechanism of action,^{4,14,15} pharmacokinetics,^{3,16} and clinical and pharmacological properties^{15,17} have been published. Cisplatin is a mutagenic and carcinogenic agent.^{18,19}

Precautions and Disclaimer

This product is for laboratory research use only. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

The product is soluble in water at 25 °C (0.253 g/100 g).¹⁵ Apparent molar solubilities are: in water (after 24 hours), 0.77 x 10⁻² in water (after 72 hours), 0.83 x 10⁻² and in 0.9% sodium chloride (after 72 hours), 0.50 x 10⁻².¹

Cisplatin is insoluble in alcohol.³ Apparent molar solubilities of cisplatin in other organic solvents are (at 25 °C):

Dimethylformamide (DMF) 6.67 x 10⁻²Dimethylacetamide (DMA) 6.00 x 10⁻²and Dimethyl sulfoxide (DMSO) 116.6 x 10⁻².¹

Cisplatin is soluble in DMF (approximately 16.6 mg/ml) with gentle heating.

The cis-trans isomerization of cisplatin occurs in aqueous solution (in the absence of catalysts) and is increased at elevated temperatures.^{15,20} The reversible substitution of water for chloride ions in cisplatin may be the primary reason for its decomposition in aqueous solutions.³

Even though cis-platin is soluble in DMSO, the use of DMSO to dissolve *cis*- or *trans*-diamminedichloro-platinum (DDP) in biological studies is **strongly** discouraged. The DMSO inserts itself into the ligand.²⁷

Stability of cisplatin in aqueous solutions was enhanced by increasing the sodium chloride (NaCl) concentration to 0.9% and was adversely affected in alkaline solutions such as sodium bicarbonate solutions.^{3,21} A solution of cisplatin in 0.9% NaCl was stable for 24 hours (after an initial loss of 3% in less than one hour) at room temperature protected from the light.^{3,22} The increased stability in NaCl solutions may be due to the availability of excess chloride ions.³ Solutions of <0.6 mg/ml in NaCl should be stored at 2-8 °C and protected from light to avoid precipitation.²² In addition to isotonic saline solution, cisplatin was stable in isotonic saline containing 5% dextrose in lactated Ringers solution.²¹ A 0.1% solution in 0.9% sodium chloride solution has a pH range of 4.5-6.0 immediately after preparation.³ It is recommended to prepare all solutions fresh and protect from light.

Cisplatin is stable in the presence of mannitol, glucose, benzyl alcohol, and parabens.^{1,21,23} Cisplatin reacts with nucleophiles such as bisulfite (an antioxidant),^{24,25} methionine, methionine containing peptides, and other compounds containing divalent sulfur.¹ Cisplatin (1 mg/ml) also reacts with aluminum producing a gas and black precipitate.²⁶

References

- Riley, C.M., and Sternson, L.A., "Cisplatin" Analytical Profiles of Drug Substances, **14**, 77, (1985).
- Berners-Price, S.J. et al., J. Chem. Soc. Chem. Commun. **789**, 1992.
- Martindale, The Extra Pharmacopoeia, 31st ed. p. 552.
- Chabner, B.A. et al. "Antineoplastic Agents" Chapter 51, p. 1233 in Goodman and Gilman's The Pharmacological Basis of Therapeutics, 9th ed. Hardman, J.G. et al. Eds. McGraw-Hill, New York, 1996.
- Rosenberg, B. et al., Nature, **222**, 385, (1969).
- Ushay, H.M. et al., Biochem., **20**, 3744, (1981).
- Peyrot, V. et al., Biochem. Pharmacol., **35**, 371, (1986).
- Boekelheide, K. et al., Toxicol. Appl. Pharmacol., **116**, 146, (1992).
- Kopf-Maier, P., and Muhlhausen, S.K., Chem. Biol. Interact., **82**, 295, (1992).
- Zaffaroni, N. et al., Br. J. Cancer, **77**, 1378, (1998).
- Henkels, K.M., and Turchi, J.J., Cancer Res., **57**, 4488, (1997).
- Speelmans, G. et al., Biochem., **36**, 10545, (1997).
- Daugaard, G. et al., J. Pharmacol. Exp. Ther., **244**, 1081, (1988).
- Lippard, S.J., Science, **218**, 1075, (1982).
- The Merck Index, 12th ed., Entry # 2378.
- Benet, L.Z. et al. "Design and Optimization of Dosage Regimens: Pharmacokinetic Data" Appendix II, 1707, in Goodman and Gilman's The Pharmacological Basis of Therapeutics, 9th ed. Hardman, J.G. et al. Eds. McGraw-Hill, New York, 1996.
- DeConti, R.C. et al., Cancer Res., **33**, 1310, (1973).
- Chibber, R. et al., Eur. J. Cancer and Clin. Oncol., **25**, 27, (1989).
- Burnouf, D. et al., Proc. Natl Acad. Sci., **87**, 6087, (1990).
- Breusova-Baidala, Y. et al., Izv. Akad. Nauk SSSR, Ser. Khim., **6**, 1239, 1974 (in Russian). (Chemical Abstracts, 78, 987089).
- Hincal, A.A. et al., J. Parenter. Drug Assoc., **33**, 107, (1979).
- Greene, R.F. et al., Am. J. Hos. Pharm., **36**, 38, (1979).
- LaFollette, J.M. et al., Am J. Hosp. Pharm., **42**, 2652, (1985).
- Hussain, A.A. et al., J. Pharm. Sci., **69**, 364, (1980).
- Garren, K.W., and Repta, A.J., Int'l J. Pharmaceutics, **24**, 91, (1985).
- Ogawa, G.S. et al., Am. J. Hosp. Pharm., **42**, 1042, (1985).
- Sundquist, W.I., et al., Solvolysis Reactions of *cis*- and *trans*-Diamminedichloroplatinum(II) in Dimethyl Sulfoxide. Structural characterization and DNA binding of *trans*-[Pt(NH₃)₂(Me₂SO)Cl]⁺. Inorg. Chem., **26**, 1524-1528 (1987).

ARO/MAM 5/07

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.