

Earle's Balanced Salts

Component	E2888 [Liquid] g/L	E3024 [Liquid] g/L	E7510 [Liquid] g/L
<u>Inorganic salts</u>			
CaCl ₂ • 2H ₂ O	0.265	0.265	2.65
MgSO ₄ (anhyd)	0.09767	0.09767	0.9767
KCl	0.4	0.4	4
NaHCO ₃	2.2	2.2	—
NaCl	6.8	6.8	68
Na ₂ HPO ₄ (anhyd)	0.122	0.122	1.22
<u>Sugars</u>			
D-Glucose	1	1	10
<u>Other</u>			
Phenol red • Na	0.011	—	0.11
<u>Add</u>			
NaHCO ₃	—	—	2.2 at 1x