


Fit to clean

New Extran® 5L canisters
for lab dishwashers




The life science business of
Merck KGaA, Darmstadt,
Germany operates as
MilliporeSigma in the
U.S. and Canada.

Supelco®
Analytical Products

Extran[®] in NEW 5L canister

Experience the new safety and simplicity perfectly aligned to your lab dishwasher

Optimal size and design



Allows direct connection to most lab dishwashers



Extran[®] – Eliminate hassles and hazards

Careful cleaning is essential in every lab. But it can pose risks and challenges, like accidental contact with harmful cleaner concentrates, frequent refilling of cleaning agents, or costly dishwasher repairs.

That's why we now also offer our **powerful, non-toxic, residue-free and biodegradable Extran[®] cleaners in new 5L canisters** – making them the **ideal fit for lab** dishwashers. Discover quality perfected for your intended use.



Explore the complete range of Extran[®] and download our "Cleaner" brochure SigmaAldrich.com/extran-br

Top label with contents & essential safety information

Order no. 1.40118.5001

Exact product designation (EN) **Extran® AP 18**
liquid, mild alkaline concentrate

Color coded with clear designation acc. to pH of the cleaner
Blue = alkaline
Red = acidic

Hazard classification / pictogram

RECEIVED: _____

OPENED: _____

EXPDATE: _____

Space for individual documentation with permanent marker

H+P statements & handling instructions

Ordering Information

Product	Content	Packaging	Ord. No.
Extran® AP 17 liquid alkaline	5 L	PE canister	1.40006.5001
Extran® AP 18 liquid mild alkaline	5 L	PE canister	1.40118.5001
Extran® AP 21 acidic with phosphoric acid	5 L	PE canister	1.07559.5001
Extran® AP 22 acidic with citric acid	5 L	PE canister	1.07561.5001

The multi-purpose cleaner with multiple benefits

Safe

- Connected canister avoids contact with harmful acidic or alkaline cleaners
- Product and safety information always visible on top label
- Top label with blank space for quick documentation of individual GLP handling dates and unique color coding acc. to pH

Convenient

- Canister fits in most dishwashers, with standard thread for direct connection
- No need to clean dishwasher canisters or refill from smaller bottles
- Easy to use thanks to light weight and large handle

Economical

- Avoid costly repairs from sediment buildup in refilled canisters
- Save resources through direct connection using concentrates

Ecological

- Extran® cleaners contain non-toxic, biodegradable active ingredients



Supelco®

Analytical Products

MilliporeSigma
290 Concord Road
Billerica, MA 01821

SigmaAldrich.com/extran-br

© 2020 Merck KGaA, Darmstadt, Germany and its affiliates. All Rights Reserved.
MilliporeSigma, the vibrant M, Supelco and Extran are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners.
Detailed information on trademarks is available via publicly accessible resources.

Lit. No. MS_FL5047EN Ver. 1.0 01/2020