

The Anatomy of GenElute™-E


single spin DNA AND RNA purification



Silica-based DNA and RNA purification, with all of its spins and washes, can be tedious.


Discover how GenElute™-E Single Spin technology enables nucleic acid purification in a single step!

GET READY...



A quick SmartLyse™ treatment prepares the DNA sample for loading onto a packed column for a unique use of our resin.

SET...



How did we simplify this process? Negative chromatography! This allows our columns to absorb and retain sample contaminants while allowing nucleic acids to flow through the column. Fewer steps and reduced waste.

GO!



In no time at all, the purified nucleic acid is ready for downstream analysis.

**Better results. Fewer steps.
Just elute. Learn more at
SigmaAldrich.com/singlespin**