

# Tumor-Associated Antigen Panel Cell Lines for Testing of Novel Immunotherapies



Revolutionizing Immuno-Oncology Research with Biologically-Relevant Cellular Assays

We have utilized CompoZr® zinc finger nuclease (ZFN) technology and MISSION® lentiviral technology to create functional cell-based potency assays that express individual tumor antigens at multiple levels in biologically-relevant cell lines. These cell lines can be used throughout the entire drug discovery process and would be ideal cell-based assays for CAR/TCR therapeutic screening, checkpoint inhibitor development, or target cell lines for potency assays in clinical development.





## Tumor Associated Antigen Panel Cell Lines

- Three gene-edited cell lines that express individual tumor antigens at KO, low, and high levels
- KO cell line generation generated with ZFN technology to enable clear IP/licensing path for potential use in clinical development
- Stable antigen expression at low/high levels generated using lentiviral technology
- Customized tumor antigen/cell line combination is available through Cell Design Studio® Services



## Tumor Associated Antigen Panel Cell Lines



### Order Information

Cat. No.	Product Name
ATG001-1KT	Raji CD19
ATG002-1KT	A549 EGFR
ATG004-1KT	A375 NY-ESO-1
A549 EGFRvIII	*A549 EGFRvIII

\*A549 EGFRvIII have KO and high expressing tumor antigens only

For more information visit:

[SigmaAldrich.com/ImmunoOncology](https://SigmaAldrich.com/ImmunoOncology)

[SigmaAldrich.com/CellDesignStudio](https://SigmaAldrich.com/CellDesignStudio)

To place an order or receive technical assistance in the U.S. and Canada, call toll-free 1-800-645-5476  
For other countries across Europe and the world, please visit: [SigmaAldrich.com/offices](https://SigmaAldrich.com/offices)  
For Technical Service, please visit: [SigmaAldrich.com/techservice](https://SigmaAldrich.com/techservice)

MilliporeSigma  
400 Summit Drive  
Burlington, MA 01803

[SigmaAldrich.com](https://SigmaAldrich.com)

© 2020 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. MilliporeSigma, the vibrant M, Mission, Compozr, Cell Design Studio and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

