

Millipore®

Preparation, Separation,
Filtration & Monitoring Products



Mobius® Breez Microbioreactor

A 2 mL Automated Perfusion
Cell Culture Platform





- Run more experiments in parallel for higher efficiency
- Reliable run-to-run performance, obtain high-quality perfusate on demand
- Freedom from hands-on work with automated control and monitoring
- Operational flexibility allows you to save on development time and costs
- Easy to set up and use, allowing you to focus on your goals

Working Volume	2 mL
# of Controllers	4 (independent)
Flow Rate	0.1 – 10 VVD
Sample Volume	50 µL (min.)
Closed Loop Control	pH, DO, T, density-OD

Do more with less and accelerate your path to process intensification

- ✓ Perfusion process development
- ✓ Cell line development
- ✓ Media screening and optimization



 **Ask your sales rep about a trial!**

MilliporeSigma is the U.S. and Canada Life Science business of Merck KGaA, Darmstadt, Germany.

© 2024 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. MilliporeSigma, the vibrant M, Millipore and Mobius are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners.

Lit. No. MS_FL13346EN V1 02/2024

**MILLIPORE
SIGMA**