



3050 Spruce Street  
Saint Louis, Missouri 63103 USA  
Telephone 800-325-5832 • (314) 771-5765  
Fax (314) 286-7828  
email: techserv@sia.com  
sigma-aldrich.com

## Product Information

### 40% UREA SOLUTION

Product Number **U7129**

#### Product Description

A filter sterilized urea solution recommended for the detection of urease activity.

#### Components (5 mls per vial)

##### Item

|                 |          |
|-----------------|----------|
| Urea            | 2.00 g   |
| Distilled Water | 5.00 mls |

#### Preparation Instructions

Warm the refrigerated 40% Urea Solution to room temperature. Aseptically add 5 mls to 95 mls of sterile molten Urea Agar Base (Christensen). Mix well and dispense in sterile tubes.

#### Storage

Store the 40% Urea Solution at 2-8°C.

#### Precautions and Disclaimer

For laboratory use only. Not for drug, household or other uses.

U7129 6/99

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.